www.neobiolab.com info@neobiolab.com 888.754.5670, +1 617.500.7103 United States 0800.088.5164, +44 020.8123.1558 United Kingdom

UNG

Description: E. Coli Uracil DNA Glycosilase (UNG) catalyses the release of free Uracil from Uracil-containing DNA. UNG efficiently hydrolyzes uracil from signle-stranded or double-stranded DNA, but not from oligomers (6 fewer bases).

Catalog #:ENPS-359

For research use only.

Synonyms: Uracil DNA Glycosilase, Uracil DNA Glycosylase, UNG.

Source: Escherichia Coli strain that carries the UNG gene from E.coli.

Physical Appearance: Sterile Filtered colorless solution.

Formulation:

UNG solution in 10mM Tris-Hcl (pH-7.4 at 25°C), 50mM KCl, 1mM DTT, 0.1mM EDTA, 0.1 mg/ml BSA and 50% glycerol.

Stability:

Uracil DNA Glycosilase although stable at 15°C for 1 week, should be stored desiccated below -18°C. Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drµgs, agricultural or pesticidal products, food additives or household chemicals.

To place an order, please Click HERE.



