



**Product Name** Fas/APO-1

**Catalog Number** MF0012 (96 Tests)

**Assay Time** 90 minutes

Sample 100 ul

Conjugate HRP

**Instrument** Microlate Reader

**Precautions** ELISA

**Product Overview** The ELISA (Enzyme-Linked Immunosorbent Assay) kit is an in vitro enzyme-

linked immunosorbent assay for the quantitative measurement of samples in

serum, plasma, cell culture supernatants and urine.

**Intended Use** This Mouse Fas/APO-1 ELISA Kit is intended for laboratory research use only

and not for use in diagnostic or therapeutic procedures. The stop solution changes color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of Mouse Fas/APO-1 in the sample, this Mouse Fas/APO-1 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of optical density versus Mouse Fas/APO-1 concentration. The concentration of the samples is then determined by comparing the O.D. of the samples to the

standard curve.

**Storage Instructions** Store at 2-8℃

**Safety Notes** 1. This kit contains 0.1% Proclin 300 as a preservative and does not contain

sodium azide. Proclin 300 is a relatively safe preservative; however, avoid ingestion and contact with eyes, skin, or mucous membranes. In case of contact, rinse the affected area with plenty of water. Observe all federal, state,

and local regulations for disposal.

2. All blood components and biological materials should be handled as potentially hazardous. Follow universal precautions as established by the Centers for Disease Control and Prevention and by the Occupational Safety and

Health Administration when handling and disposing of infectious agents.

FOR RESEARCH USE ONLY. NOT FOR THERAPEUTIC OR DIAGNOSTIC APPLICATIONS. READ THROUGH ALL PROCEDURES BEFORE USE.

Tel: +1(888) 733-6849 (Toll Free) or +1(617) 299-7367 Fax: +1(888) 733-6849

Email: sales@NeoScientific.com https://www.neobiolab.com