

YAP1

Reactivity: Human

Tested applications: WB IHC ICC IF IP FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:50 - 1:200 IF 1:50 - 1:100
IP 1:20 - 1:50 FC 1:20 - 1:50

Calculated MW: 65 kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human YAP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

YAP; YKI; COB1; YAP2; YAP65;

Catalog #: A11326

Antibody Type:

Monoclonal Antibody

Species: Rabbit

Gene ID: 10413

Isotype: IgG

Swiss Prot: P46937

Purity: Affinity purification

For research use only.

Background:

This gene encodes a downstream nuclear effector of the Hippo signaling pathway which is involved in development, growth, repair, and homeostasis. This gene is known to play a role in the development and progression of multiple cancers as a transcriptional regulator of this signaling pathway and may function as a potential target for cancer treatment. Alternative splicing results in multiple transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).