

YAP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 54kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human YAP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AP; YAP2; YAP65; YKI

Catalog #: A11430

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10413

Isotype: IgG

Swiss Prot: P46937

Purity: Affinity purification

For research use only.

Background:

This gene encodes a downstream nuclear effector of the Hippo signaling pathway which is involved in development, growth, repair, and homeostasis. This gene is known to play a role in the development and progression of multiple cancers as a transcriptional regulator of this signaling pathway and may function as a potential target for cancer treatment. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2013]

To place an order, please [Click HERE](#).