

PRDM14

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP ChIP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200 IP 1:20 - 1:100

ChIP 1:20 - 1:100

Calculated MW: 64kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PRDM14

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

PFM11;

Catalog #: A5543

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 63978

Isotype: IgG

Swiss Prot: Q9GZV8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the PRDI-BF1 and RIZ homology domain containing (PRDM) family of transcriptional regulators. The encoded protein may possess histone methyltransferase activity and plays a critical role in cell pluripotency by suppressing the expression of differentiation marker genes. Expression of this gene may play a role in breast cancer.

To place an order, please [Click HERE](#).