

L3MBTL1

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 84kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human L3MBTL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

L3MBTL; ZC2HC3; H-L(3)MBT; dJ138B7.3;

Catalog #: A7090

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26013

Isotype: IgG

Swiss Prot: Q9Y468

Purity: Affinity purification

For research use only.

Background:

This gene represents a polycomb group gene. The encoded protein functions to regulate gene activity, likely via chromatin modification. The encoded protein may also be necessary for mitosis. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).