

BMP15

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 45kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BMP15

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ODG2; POF4; GDF9B;

Catalog #: A7321

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9210

Isotype: IgG

Swiss Prot: O95972

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the bone morphogenetic protein family which is part of the transforming growth factor-beta superfamily. The transforming growth factor-beta superfamily includes large families of growth and differentiation factors. It is thought that this protein may be involved in oocyte maturation and follicular development as a homodimer or by forming heterodimers with a related protein, Gdf9. Defects in this gene are the cause of ovarian dysgenesis 2.

To place an order, please [Click HERE](#).