

LHX6

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 40kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LHX6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LHX6.1;

Catalog #: A7364

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26468

Isotype: IgG

Swiss Prot: Q9UPM6

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of a large protein family that contains the LIM domain, a unique cysteine-rich zinc-binding domain. The encoded protein may function as a transcriptional regulator and may be involved in the control of differentiation and development of neural and lymphoid cells. Two alternatively spliced transcript variants encoding distinct isoforms have been described for this gene. Alternatively spliced transcript variants have been identified, but their biological validity has not been determined.

To place an order, please [Click HERE](#).