

BLK

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 58kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BLK

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MODY11;

Catalog #: A7427

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 640

Isotype: IgG

Swiss Prot: P51451

Purity: Affinity purification

For research use only.

Background:

This gene encodes a nonreceptor tyrosine-kinase of the src family of proto-oncogenes that are typically involved in cell proliferation and differentiation. The protein has a role in B-cell receptor signaling and B-cell development. The protein also stimulates insulin synthesis and secretion in response to glucose and enhances the expression of several pancreatic beta-cell transcription factors.

To place an order, please [Click HERE](#).