

GSDMB

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GSDMB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GSDML; PP4052; PRO2521;

Catalog #: A7474

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55876

Isotype: IgG

Swiss Prot: Q8TAX9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the gasdermin-domain containing protein family. Other gasdermin-family genes are implicated in the regulation of apoptosis in epithelial cells, and are linked to cancer. Multiple transcript variants encoding different isoforms have been found for this gene. Additional variants have been described, but they are candidates for nonsense-mediated mRNA decay (NMD) and are unlikely to be protein-coding.

To place an order, please [Click HERE](#).