

GABRA3

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 55kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human GABRA3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A7522

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2556

Isotype: IgG

Swiss Prot: P20236

Purity: Affinity purification

For research use only.

Background:

GABA is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA-A receptors, which are ligand-gated chloride channels. Chloride conductance of these channels can be modulated by agents such as benzodiazepines that bind to the GABA-A receptor. At least 16 distinct subunits of GABA-A receptors have been identified.

To place an order, please [Click HERE](#).