

VANGL1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 60kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human VANGL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LPP2; STB2; STBM2; KITENIN;

Catalog #: A7606

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 81839

Isotype: IgG

Swiss Prot: Q8TAA9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the tetraspanin family. The encoded protein may be involved in mediating intestinal trefoil factor induced wound healing in the intestinal mucosa. Mutations in this gene are associated with neural tube defects. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).