

KIF6

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 93kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human KIF6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C6orf102; dJ137F1.4; dJ188D3.1; dJ1043E3.1;

To place an order, please [Click HERE](#).

Catalog #: A7614

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 221458

Isotype: IgG

Swiss Prot: Q6ZMV9

Purity: Affinity purification

For research use only.