

## GNAI2

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 40kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human GNAI2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

GIP; GNAI2B; H\_LUCA15.1; H\_LUCA16.1;

**Catalog #:** A7676

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 2771

**Isotype:** IgG

**Swiss Prot:** P04899

**Purity:** Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is an alpha subunit of guanine nucleotide binding proteins (G proteins). The encoded protein contains the guanine nucleotide binding site and is involved in the hormonal regulation of adenylate cyclase. Several transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*