

## TFG

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 43kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human TFG

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

TF6; HMSNP; SPG57; TRKT3;

**Catalog #:** A7769

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 10342

**Isotype:** IgG

**Swiss Prot:** Q92734

**Purity:** Affinity purification

For research use only.

**Background:**

There are several documented fusion oncoproteins encoded partially by this gene. This gene also participates in several oncogenic rearrangements resulting in anaplastic lymphoma and mixoid chondrosarcoma, and may play a role in the NF-kappaB pathway. Multiple transcript variants have been found for this gene.

*To place an order, please [Click HERE](#).*