

## KCNS3

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 56kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human KCNS3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

KV9.3;

**Catalog #:** A7906

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 3790

**Isotype:** IgG

**Swiss Prot:** Q9BQ31

**Purity:** Affinity purification

For research use only.

**Background:**

Voltage-gated potassium channels form the largest and most diversified class of ion channels and are present in both excitable and nonexcitable cells. Their main functions are associated with the regulation of the resting membrane potential and the control of the shape and frequency of action potentials. The alpha subunits are of 2 types: those that are functional by themselves and those that are electrically silent but capable of modulating the activity of specific functional alpha subunits. The protein encoded by this gene is not functional by itself but can form heteromultimers with member 1 and with member 2 (and possibly other members) of the Shab-related subfamily of potassium voltage-gated channel proteins. This gene belongs to the S subfamily of the potassium channel family. Alternatively spliced transcript variants encoding the same protein have been found for this gene.

**To place an order, please [Click HERE](#).**