

## HDAC4

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 197kDa

**Observed MW:** Refer to figures

**Immunogen:**

A synthetic peptide of human HDAC4

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HDAC4; HA6116; HD4; HDAC-A; HDACA; KIAA0288 ;

*To place an order, please [Click HERE](#).*

**Catalog #:** A7951

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 9759

**Isotype:** IgG

**Swiss Prot:** P56524

**Purity:** Affinity purification

For research use only.