

USP25

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 122kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human USP25

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

USP21;

Catalog #: A7975

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29761

Isotype: IgG

Swiss Prot: Q9UHP3

Purity: Affinity purification

For research use only.

Background:

Ubiquitin (MIM 191339) is a highly conserved 76-amino acid protein involved in regulation of intracellular protein breakdown, cell cycle regulation, and stress response. Ubiquitin is released from degraded proteins by disassembly of the polyubiquitin chains, which is mediated by ubiquitin-specific proteases (USPs), such as USP25 (Valero et al., 1999)

To place an order, please [Click HERE](#).