

## USP26

---

**Reactivity:** Human

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:** 104kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human USP26

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:** A7999

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 83844

**Isotype:** IgG

**Swiss Prot:** Q9BXU7

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a member of the ubiquitin-specific processing (UBP) family of proteases and is a deubiquitinating enzyme (DUB) with His and Cys domains. It is specifically expressed in testis tissue. Mutations in this gene have been associated with Sertoli cell-only syndrome and male infertility.

*To place an order, please [Click HERE](#).*