

MLF1

Reactivity: Human

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW: 31kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MLF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A8012

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4291

Isotype: IgG

Swiss Prot: P58340

Purity: Affinity purification

For research use only.

Background:

This gene encodes an oncoprotein which is thought to play a role in the phenotypic determination of hemopoietic cells. Translocations between this gene and nucleophosmin have been associated with myelodysplastic syndrome and acute myeloid leukemia. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).