

MTCH1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 42kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MTCH1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PSAP; PIG60; CGI-64; SLC25A49;

Catalog #: A8063

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23787

Isotype: IgG

Swiss Prot: Q9NZJ7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the mitochondrial carrier family. The encoded protein is localized to the mitochondrion inner membrane and induces apoptosis independent of the proapoptotic proteins Bax and Bak. Pseudogenes on chromosomes 6 and 11 have been identified for this gene. Alternatively spliced transcript variants encoding multiple isoforms have been observed.

To place an order, please [Click HERE](#).