

## NKX3-2

---

**Reactivity:**Human Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:200 - 1:1000

**Calculated MW:**35kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human NKX3-2

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

NK3 homeobox 2; BAPX1, Homeobox protein Nkx 3.2, NK3 homeobox 2, NKX3 2, NKX3.2, NKX3-2, NKX3B;

**Background:**

This gene encodes a member of the NK family of homeobox-containing proteins. The encoded protein may play a role in skeletal development.

*To place an order, please [Click HERE](#).*

**Catalog #:**A8093

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**579

**Isotype:**IgG

**Swiss Prot:**P78367

**Purity:**Affinity purification

For research use only.