

ACRV1

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:28kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ACRV1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SP-10; SPACA2; D11S4365;

Catalog #:A8098

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:56

Isotype:IgG

Swiss Prot:P26436

Purity:Affinity purification

For research use only.

Background:

This gene encodes a testis-specific, differentiation antigen, acrosomal vesicle protein 1, that arises within the acrosomal vesicle during spermatogenesis, and is associated with the acrosomal membranes and matrix of mature sperm. The acrosomal vesicle protein 1 may be involved in sperm-zona binding or penetration. Alternatively spliced transcript variants have been described.

To place an order, please [Click HERE](#).