

DLX2

Reactivity:Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:34kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DLX2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TES1; TES-1;

Catalog #:A8410

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1746

Isotype:IgG

Swiss Prot:Q07687

Purity:Affinity purification

For research use only.

Background:

Many vertebrate homeo box-containing genes have been identified on the basis of their sequence similarity with Drosophila developmental genes. Members of the Dlx gene family contain a homeobox that is related to that of Distal-less (Dll), a gene expressed in the head and limbs of the developing fruit fly. The Distal-less (Dlx) family of genes comprises at least 6 different members, DLX1-DLX6. The DLX proteins are postulated to play a role in forebrain and craniofacial development. This gene is located in a tail-to-tail configuration with another member of the gene family on the long arm of chromosome 2.

To place an order, please [Click HERE](#).