

## GNRH2

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:100

**Calculated MW:** 13kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human GNRH2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

GnRH-II; LH-RHII;

**Catalog #:** A8424

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 2797

**Isotype:** IgG

**Swiss Prot:** O43555

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a secreted peptide hormone and member of the gonadotropin-releasing hormone (GnRH) family of proteins. The encoded protein regulates reproductive function by stimulating the production and release of the gonadotropins follicle-stimulating hormone (FSH) and luteinizing hormone (LH). The encoded protein may inhibit endometrial, ovarian, prostate, and breast cancer cell proliferation. Alternative splicing results in multiple transcript variants.

*To place an order, please [Click HERE](#).*