

NFKBIE

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:53kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human NFKBIE

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

kp

Synonym:

IKBE;

Catalog #:A8443

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4794

Isotype:IgG

Swiss Prot:O00221

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene binds to components of NF-kappa-B, trapping the complex in the cytoplasm and preventing it from activating genes in the nucleus. Phosphorylation of the encoded protein targets it for destruction by the ubiquitin pathway, which activates NF-kappa-B by making it available to translocate to the nucleus.

To place an order, please [Click HERE](#).