

CPQ

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:52kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CPQ

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LDP; PGCP;

Catalog #:A8478

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10404

Isotype:IgG

Swiss Prot:Q9Y646

Purity:Affinity purification

For research use only.

Background:

This gene encodes a metallopeptidase that belongs to the peptidase M28 family. The encoded protein may catalyze the cleavage of dipeptides with unsubstituted terminals into amino acids.

To place an order, please [Click HERE](#).