

EXO5

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 42kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human EXO5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DEM1; Exo V; hExo5; C1orf176;

To place an order, please [Click HERE](#).

Catalog #: A8507

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64789

Isotype: IgG

Swiss Prot: Q9H790

Purity: Affinity purification

For research use only.