

ZIC3

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 51kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ZIC3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HTX; HTX1; ZNF203; VACTERLX

Catalog #: A9476

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7547

Isotype: IgG

Swiss Prot: O60481

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the ZIC family of C2H2-type zinc finger proteins. This nuclear protein probably functions as a transcription factor in early stages of left-right body axis formation. Mutations in this gene cause X-linked visceral heterotaxy, which includes congenital heart disease and left-right axis defects in organs.

To place an order, please [Click HERE](#).