

TRIP12

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:220kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TRIP12

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A9958

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9320

Isotype:IgG

Swiss Prot:Q14669

Purity:Affinity purification

For research use only.

To place an order, please [Click HERE](#).