

## Phospho-PTPN11-Y580

**Reactivity:** Human Mouse Rat

**Tested applications:** IHC

**Recommended Dilution:** IHC 1:50 - 1:100

**Calculated MW:** 72kDa

**Observed MW:** Refer to Figures

**Immunogen:**

A phospho specific peptide corresponding to residues surrounding Y580 of human PTPN11

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

o

**Synonym:**

CFC; NS1; SHP2; BPTP3; PTP2C; PTP-1D; SH-PTP2; SH-PTP3;

**Catalog #:** AP0439

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 5781

**Isotype:** IgG

**Swiss Prot:** Q06124

**Purity:** Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP contains two tandem Src homology-2 domains, which function as phospho-tyrosine binding domains and mediate the interaction of this PTP with its substrates. This PTP is widely expressed in most tissues and plays a regulatory role in various cell signaling events that are important for a diversity of cell functions, such as mitogenic activation, metabolic control, transcription regulation, and cell migration. Mutations in this gene are a cause of Noonan syndrome as well as acute myeloid leukemia. Two transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*