

IFN a 2b Human, HEK

Description: IFN-a 2b Human Recombinant produced in HEK cells is a glycosylated monomer, having a total molecular weight of 16kDa. The IFN-a 2b is purified by proprietary chromatographic techniques.

Catalog #: CYPs-099

For research use only.

Synonyms: Interferon alpha 2b, IFNA, INFA2, MGC125764, MGC125765.

Source: HEK.

Physical Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Purity: Greater than 95% as observed by SDS-PAGE.

Formulation:

The IFN-a 2b was lyophilized from 1mg/ml in 1xPBS.

Stability:

Lyophilized IFN-a 2b although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution IFN-a 2b should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolabs products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

It is recommended to reconstitute the lyophilized IFN-a 2b in sterile water not less than 100

Introduction:

IFN-alpha is produced by macrophages and has antiviral activities. Interferon stimulates the production of two enzymes: protein kinase and an oligoadenylate synthetase.

Biological Activity:

The specific activity was determined by the dose-dependent cytotoxicity of the human TF-1 cell line (human erythroleukemic indicator cell line) and is typically 0.02-0.08ng/ml.

To place an order, please [Click HERE](#).