

AP3S1 Human

Description: AP3S1 Human Recombinant produced in E. coli is a single polypeptide chain containing 213 amino acids (1-193) and having a molecular mass of 23.8 kDa. AP3S1 is fused to 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: PRPS-1287

For research use only.

Synonyms: CLAPS3, Sigma3A, AP-3 complex subunit sigma-1, AP-3 complex subunit sigma-3A, Adapter-related protein complex 3 sigma-1 subunit, Clathrin-associated/assembly/adaptor protein, small 3, Sigma-3A-adaptin, Sigma-adaptin 3a.

Source: E.coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHH SSGLVPRGSH MIKAILIFNN HGKPRLSKFY
QPYSEDYQQ IIRETFHLVS RDNVCNLF EGGLIGGSD NKLIYRHYAT LYFVFCVDSS
ESELGILDLI QVFVETLDKC ENVCELDLI FHVDKVHNL AEMVMGGMVL ETNMNEIVTQ
IDAKNKLEKS EAGLAGAPAR AVSAVKNMNL PEIPRNINIG DISIKVPNLP SFK.

Purity: Greater than 90% as determined by SDS-PAGE.

Formulation:

The AP3S1 solution (1mg/1ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 10% glycerol and 1mM DTT.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

AP-3 complex subunit sigma-1 (AP3S1) belongs to the AP-3 complex, an adapter-related complex which is not clathrin-associated. AP3S1 is linked with the Golgi region in addition to more peripheral structures. AP3S1 facilitates the budding of vesicles from the Golgi membrane and involved in trafficking to lysosomes.

To place an order, please [Click HERE](#).