

DNAL4 Human

Description: DNAL4 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 125 amino acids (1-105 a.a) and having a molecular mass of 14.1kDa. DNAL4 is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: PRPS-234

For research use only.

Synonyms: Dynein light chain 4 axonemal, DNAL4, PIG27.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MGETEGKKDE ADYKRLQTFP
LVRHSDMPEE MRVETMELCV TACEKFSNNN ESAAKMIKET MDKKFGSSWH VVIGEGFGFE
ITHEVKNNLY LYFGGTLAVC VWKCS.

Purity: Greater than 95% as determined by SDS-PAGE.

Formulation:

DNAL4 protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH8.0), 20% glycerol, 0.1M NaCl and 1mM DTT.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolabs products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

DNAL4 is a member of the dynein light chain family. DNAL4 is a component of the dynein motor complex in cells that converts the chemical energy contained in ATP into the mechanical energy of movement. Dynein transports various cellular cargos by moving the length of cytoskeletal microtubules towards the minus-end of the microtubule, which is typically oriented towards the cell center.

To place an order, please [Click HERE](#).