

FYN

Reactivity: Human Mouse Rat

Tested applications: WB ICC

Recommended Dilution: WB 1:200 - 1:1000 ICC 1:50 - 1:100

Calculated MW: 59kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human FYN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FYN;SLK;SYN;p59-FYN;Tyrosine-protein kinase Fyn antibody;p59-Fyn antibody;

Catalog #: A0086

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 2534

Isotype: IgG

Swiss Prot: P06241

Purity: Affinity purification

For research use only.

Background:

Fyn is a member of the protein-tyrosine kinase oncogene family. It is a membrane-associated tyrosine kinase that has been implicated in the control of cell growth. The protein associates with the p85 subunit of phosphatidylinositol 3-kinase and interacts with the fyn-binding protein.

To place an order, please [Click HERE](#).