

IL2

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 18kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human IL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

IL2; IL-2; TCGF; lymphokine ;

Catalog #: A0302

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3558

Isotype: IgG

Swiss Prot: P60568

Purity: Affinity purification

For research use only.

Background:

IL-2 is a T cell stimulatory cytokine best known for inducing T cell proliferation and NK cell proliferation and activation (1,2). IL-2 also promotes peripheral development of Treg cells (4, 5). Conversely, IL-2 is involved in activation induced cell death (AICD) that is observed post T cell expansion by increasing levels of Fas on CD4+ T cells (3). The effects of IL-2 are mediated through a trimeric receptor complex consisting of IL-2R, IL-2R and the common gamma chain, c (1,2). IL-2 R binds exclusively to IL-2 with low affinity and increases binding affinity of the whole receptor complex including IL-2R and c subunits. IL-2R is also used by IL-15 (1,2). The common c is used by other cytokines including IL-4, IL-7, IL-9, IL-15, and IL-21 (1,2). Binding of IL-2 initiates signaling cascades involving Jak1, Jak3, Stat5 and the PI3K/Akt pathways (1,2).

To place an order, please [Click HERE](#).