

LMNA

Reactivity:Human

Tested applications:WB

Calculated MW:664

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human LMNA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A0464

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4000

Isotype:IgG

Swiss Prot:P02545

Purity:Affinity purification

For research use only.