

IL6ST

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW:104kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human IL6ST

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD130; GP130; CDW130; IL-6RB;

Catalog #:A0469

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3572

Isotype:IgG

Swiss Prot:P40189

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a signal transducer shared by many cytokines, including interleukin 6 (IL6), ciliary neurotrophic factor (CNTF), leukemia inhibitory factor (LIF), and oncostatin M (OSM). This protein functions as a part of the cytokine receptor complex. The activation of this protein is dependent upon the binding of cytokines to their receptors. vIL6, a protein related to IL6 and encoded by the Kaposi sarcoma-associated herpesvirus, can bypass the interleukin 6 receptor (IL6R) and directly activate this protein. Knockout studies in mice suggest that this gene plays a critical role in regulating myocyte apoptosis. Alternatively spliced transcript variants encoding distinct isoforms have been described. A related pseudogene has been identified on chromosome 17.

To place an order, please [Click HERE](#).