

TUBB

Reactivity: Human Mouse Rat

Tested applications: WB IF

Recommended Dilution: WB 1:2000 - 1:10000 IF 1:50 - 1:100

Calculated MW: 50kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TUBB

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A0481

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 203068

Isotype: IgG

Swiss Prot: P07437

Purity: Affinity purification

For research use only.