

## SKP1

**Reactivity:** Human

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 19kDa

**Observed MW:** Refer to Figures

**Immunogen:**

A synthetic peptide of human SKP1

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

d

**Synonym:**

EMC19; MGC34403; OCP-II; OCP2; SKP1A; TCEB1L; p19A

**Catalog #:** A0554

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 6500

**Isotype:** IgG

**Swiss Prot:** P63208

**Purity:** Affinity purification

For research use only.

**Background:**

SKP1(S-phase kinase-associated protein 1) is also named as EMC19, OCP2, SKP1A, TCEB1L, belongs to the SKP1 family. It is an essential component of the SCF (SKP1-CUL1-F-box protein) ubiquitin ligase complex, which mediates the ubiquitination of proteins involved in cell cycle progression, signal transduction and transcription. This antibody is specific to SKP1.

*To place an order, please [Click HERE](#).*