

ITK

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 72kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human ITK

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EMT; LYK; MGC126257; MGC126258; PSCTK2

Catalog #: A0566

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3702

Isotype: IgG

Swiss Prot: Q08881

Purity: Affinity purification

For research use only.

Background:

This gene encodes an intracellular tyrosine kinase expressed in T-cells. The protein contains both SH2 and SH3 domains which are often found in intracellular kinases. It is thought to play a role in T-cell proliferation and differentiation. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).