

EDN1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 24kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human EDN1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

a

Synonym:

EDN1;ET1;HDLCQ7 ;

Catalog #: A0686

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1906

Isotype: IgG

Swiss Prot: P05305

Purity: Affinity purification

For research use only.

Background:

EDN1, also named as Endothelin-1, PPET1 and ET-1, belongs to the endothelin/sarafotoxin family. Endothelins are endothelium-derived vasoconstrictor peptides. Emerging basic science, and animal and human data all suggest that EDN1 is a potentially important contributor in the pathobiology of fibrosing disorders, including those that affect the lung (PMID:20055532).

To place an order, please [Click HERE](#).