

ANGPTL3

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:54kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human ANGPTL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

a

Synonym:

ANGPTL3;ANGPT5 ;

Catalog #:A0689

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:27329

Isotype:IgG

Swiss Prot:Q9Y5C1

Purity:Affinity purification

For research use only.

Background:

Angiotensin-like protein 3 (Angptl3) functions as a potent lipoprotein lipase inhibitor and is an important component of plasma triglyceride homeostasis. Angptl3 also plays a role in adipose formation and angiogenesis through its interaction with integrin

To place an order, please [Click HERE](#).