

TSC1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 130kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human TSC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

LAM; TSC; TSC1; KIAA0243; Hamartin; Tuberous sclerosis 1 protein;

Catalog #: A0720

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7248

Isotype: IgG

Swiss Prot: Q92574

Purity: Affinity purification

For research use only.

Background:

This gene encodes a growth inhibitory protein thought to play a role in the stabilization of tuberlin. Mutations in this gene have been associated with tuberous sclerosis. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).