

JAK3

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:115kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human JAK3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

JAKL; LJAK; JAK-3; L-JAK; JAK3_HUMAN;

Background:

The protein encoded by this gene is a member of the Janus kinase (JAK) family of tyrosine kinases involved in cytokine receptor-mediated intracellular signal transduction. It is predominantly expressed in immune cells and transduces a signal in response to its activation via tyrosine phosphorylation by interleukin receptors. Mutations in this gene are associated with autosomal SCID (severe combined immunodeficiency disease).

To place an order, please [Click HERE](#).

Catalog #:A0748

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3718

Isotype:IgG

Swiss Prot:P52333

Purity:Affinity purification

For research use only.