

RAN

Reactivity: Human Mouse

Tested applications: WB IHC FC

Calculated MW: 216

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human RAN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A0819

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5901

Isotype: IgG

Swiss Prot: P62826

Purity: Affinity purification

For research use only.