

## PRKAB1

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC ICC IP

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 ICC 1:20 - 1:50 IP 1:50 - 1:200

**Calculated MW:**38kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human PRKAB1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

c

**Synonym:**

AMPK; HAMPKb;

**Catalog #:**A0846

**Antibody Type:**

Monoclonal Antibody

**Species:**Mouse

**Gene ID:**5564

**Isotype:**IgG

**Swiss Prot:**Q9Y478

**Purity:**Affinity purification

For research use only.

**Background:**

Moesin (for membrane-organizing extension spike protein) is a member of the ERM family which includes ezrin and radixin. ERM proteins appear to function as cross-linkers between plasma membranes and actin-based cytoskeletons. Moesin is localized to filopodia and other membranous protrusions that are important for cell-cell recognition and signaling and for cell movement.

*To place an order, please [Click HERE](#).*