

BCAR1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW:93kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human BCAR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bq

Synonym:

CAS; CAS1; CASS1; CRKAS; P130Cas

Catalog #:A0855

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9564

Isotype:IgG

Swiss Prot:P56945

Purity:Affinity purification

For research use only.

Background:

BCAR1, or CAS, is an Src (MIM 190090) family kinase substrate involved in various cellular events, including migration, survival, transformation, and invasion (Sawada et al., 2006 [PubMed 17129785]).[supplied by OMIM, May 2009]

To place an order, please [Click HERE](#).