

BCAR1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW: 93kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human BCAR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bq

Synonym:

CAS; CAS1; CASS1; CRKAS; P130Cas

Catalog #: A0855

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9564

Isotype: IgG

Swiss Prot: P56945

Purity: Affinity purification

For research use only.

Background:

BCAR1, or CAS, is an Src (MIM 190090) family kinase substrate involved in various cellular events, including migration, survival, transformation, and invasion (Sawada et al., 2006 [PubMed 17129785]).[supplied by OMIM, May 2009]

To place an order, please [Click HERE](#).