

YY1

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:68kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human YY1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DELTA; NF-E1; UCRBP; INO80S; YIN-YANG-1;

Catalog #:A0856

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7528

Isotype:IgG

Swiss Prot:P25490

Purity:Affinity purification

For research use only.

Background:

YY1 is a ubiquitously distributed transcription factor belonging to the GLI-Kruppel class of zinc finger proteins. The protein is involved in repressing and activating a diverse number of promoters. YY1 may direct histone deacetylases and histone acetyltransferases to a promoter in order to activate or repress the promoter, thus implicating histone modification in the function of YY1.

To place an order, please [Click HERE](#).